




Autoanalyzer Tubing - Standard Configurations and Formulations

Two-Stop Configuration				Collar Spacing	Tygon® R-3603	Tygon® F-4040-A	Fluran®	Pharmed® BPT	Tygon® 3350
				Length (A)	15 in. (381 mm)	15 in. (381 mm)	7 in. (178 mm)	15 in. (381 mm)	15 in. (381 mm)
				Collar Distance (B)	5.98 in. (152 mm)	5.51 in. (140 mm)	5.00 in. (140 mm)	5.98 in. (152 mm)	5.00 in. (127 mm)
ID inches	ID mm	Wall Thickness inches	Wall Thickness mm	Color Code	Tygon® R-3603	Tygon® F-4040-A	Fluran®	Pharmed® BPT	Tygon® 3350
0.005	0.13	0.035	0.90	ORANGE / BLACK	A053400				
0.008	0.19	0.035	0.90	ORANGE / RED	A053401				
0.010	0.25	0.035	0.90	ORANGE / BLUE	A053402	A053303		A053903	
0.015	0.38	0.035	0.90	ORANGE / GREEN	A053403	A053304		A053904	
0.017	0.44	0.035	0.90	GREEN / YELLOW	A053404				
0.020	0.51	0.035	0.90	ORANGE / YELLOW	A053405	A053305	A065105	A053905	
0.022	0.57	0.035	0.90	WHITE / YELLOW	A053406				
0.025	0.64	0.035	0.90	ORANGE / WHITE	A053407	A053306	A065106	A053906	A049707
0.030	0.76	0.033	0.85	BLACK / BLACK	A053408	A053307	A065107	A053907	A049708
0.035	0.89	0.033	0.85	ORANGE / ORANGE	A053409	A053308	A065108	A053908	A049709
0.037	0.95	0.033	0.85	WHITE / BLACK	A053410				
0.040	1.02	0.033	0.85	WHITE / WHITE	A053411	A053309	A065109	A053909	A049710
0.043	1.09	0.033	0.85	WHITE / RED	A053412				
0.045	1.14	0.033	0.85	RED / RED	A053413	A053310	A065110	A053910	A049711
0.048	1.22	0.033	0.85	RED / GRAY	A053414				
0.051	1.30	0.033	0.85	GRAY / GRAY	A053415	A053311	A065111	A053911	A049712
0.056	1.42	0.033	0.85	YELLOW / YELLOW	A053416	A053312	A065112	A053912	A049713
0.060	1.52	0.033	0.85	YELLOW / BLUE	A053417	A053319	A065119	A053913	A049714
0.065	1.65	0.033	0.85	BLUE / BLUE	A053418	A053313	A065113	A053919	A049715
0.069	1.75	0.033	0.85	BLUE / GREEN	A053419				
0.073	1.85	0.033	0.85	GREEN / GREEN	A053420	A053314	A065114	A053914	A049716
0.081	2.06	0.033	0.85	PURPLE / PURPLE	A053421	A053315	A065115	A053915	A049717
0.090	2.29	0.033	0.85	PURPLE / BLACK	A053422	A053316	A065116	A053916	A049718
0.100	2.54	0.033	0.85	PURPLE / ORANGE	A053423	A053317	A065117	A053917	A049719
0.110	2.79	0.033	0.85	PURPLE / WHITE	A053424	A053318	A065118	A053918	A049720
0.125	3.17	0.033	0.85	BLACK / WHITE	A053425				

Three-Stop Configuration				Collar Spacing	Tygon® R-3603	Tygon® F-4040-A	Fluran®	Pharmed® BPT	Tygon® 3350
				Length (A)	15 in. (381 mm)	15 in. (381 mm)	15 in. (381 mm)	15 in. (381 mm)	15 in. (381 mm)
				Collar Distance (B)	5.98 in. (152 mm)	6.69 in. (170 mm)	6.10 in. (155 mm)	5.98 in. (152 mm)	6.50 in. (165 mm)
ID inches	ID mm	Wall Thickness inches	Wall Thickness mm	Color Code	Tygon® R-3603	Tygon® F-4040-A	Fluran®	Pharmed® BPT	Tygon® 3350
0.005	0.13	0.035	0.90	ORANGE / BLACK	A053500				
0.008	0.19	0.035	0.90	ORANGE / RED	A053501				
0.010	0.25	0.035	0.90	ORANGE / BLUE	A053502	A054303		A054003	
0.015	0.38	0.035	0.90	ORANGE / GREEN	A053503	A054304		A054004	
0.017	0.44	0.035	0.90	GREEN / YELLOW	A053504				
0.020	0.51	0.035	0.90	ORANGE / YELLOW	A053505	A054305	A065205	A054005	
0.022	0.57	0.035	0.90	WHITE / YELLOW	A053506				
0.025	0.64	0.035	0.90	ORANGE / WHITE	A053507	A054306	A065206	A054006	A049807
0.030	0.76	0.033	0.85	BLACK / BLACK	A053508	A054307	A065207	A054007	A049808
0.035	0.89	0.033	0.85	ORANGE / ORANGE	A053509	A054308	A065208	A054008	A049809
0.037	0.95	0.033	0.85	WHITE / BLACK	A053510				
0.040	1.02	0.033	0.85	WHITE / WHITE	A053511	A054309	A065209	A054009	A049810
0.043	1.09	0.033	0.85	WHITE / RED	A053512				
0.045	1.14	0.033	0.85	RED / RED	A053513	A054310	A065210	A054010	A049811
0.048	1.22	0.033	0.85	RED / GRAY	A053514				
0.051	1.30	0.033	0.85	GRAY / GRAY	A053515	A054311	A065211	A054011	A049812
0.056	1.42	0.033	0.85	YELLOW / YELLOW	A053516	A054312	A065212	A054012	A049813
0.060	1.52	0.033	0.85	YELLOW / BLUE	A053517	A054319	A065219	A054019	A049814
0.065	1.65	0.033	0.85	BLUE / BLUE	A053518	A054313	A065213	A054013	A049815
0.069	1.75	0.033	0.85	BLUE / GREEN	A053519				
0.073	1.85	0.033	0.85	GREEN / GREEN	A053520	A054314	A065214	A054014	A049816
0.081	2.06	0.033	0.85	PURPLE / PURPLE	A053521	A054315	A065215	A054015	A049817
0.090	2.29	0.033	0.85	PURPLE / BLACK	A053522	A054316	A065216	A054016	A049818
0.100	2.54	0.033	0.85	PURPLE / ORANGE	A053523	A054317	A065217	A054017	A049819
0.110	2.79	0.033	0.85	PURPLE / WHITE	A053524	A054318	A025218	A054018	A049820
0.125	3.17	0.033	0.85	BLACK / WHITE	A053525				

Nominal Tubing Without Collars				Tygon® R-3603	Tygon® F-4040-A	Fluran®	Pharmed® BPT	Tygon® 3350	
				Length (A)	33 ft. (10 m)	33 ft. (10 m)	33 ft. (10 m)	10 ft. (3m)	33 ft. (3 m)
ID		Wall Thickness		Tygon® R-3603	Tygon® F-4040-A	Fluran®	Pharmed® BPT	Tygon® 3350	
inches	mm	inches	mm						
0.005	0.13	0.035	0.90	A053400L					
0.008	0.19	0.035	0.90	A053401L					
0.010	0.25	0.035	0.90	A053402L	A053303F		A053903F		
0.015	0.38	0.035	0.90	A053403L	A053304F		A053904F		
0.017	0.44	0.035	0.90	A053404L					
0.020	0.51	0.035	0.90	A053405L	A053305F	A053805	A053905F	A049706F	
0.022	0.57	0.035	0.90	A053406L					
0.025	0.64	0.035	0.90	A053407L	A053306F	A053806	A053906F	A049707F	
0.030	0.76	0.033	0.85	A053408L	A053307F	A053807	A053907F	A049708F	
0.035	0.89	0.033	0.85	A053409L	A053308F	A053808	A053908F	A049709F	
0.037	0.95	0.033	0.85	A053410L					
0.040	1.02	0.033	0.85	A053411L	A053309F	A053809	A053909F	A049710F	
0.043	1.09	0.033	0.85	A053412L					
0.045	1.14	0.033	0.85	A053413L	A053310F	A053810	A053910F	A049711F	
0.048	1.22	0.033	0.85	A053414L					
0.051	1.30	0.033	0.85	A053415L	A053311F	A053811	A053911F	A049712F	
0.056	1.42	0.033	0.85	A053416L	A053312F	A053812	A053912F	A049713F	
0.060	1.52	0.033	0.85	A053417L	A053319F	A053819	A053919F	A049714F	
0.065	1.65	0.033	0.85	A053418L	A053313F	A053813	A053913F	A049715F	
0.069	1.75	0.033	0.85	A053419L					
0.073	1.85	0.033	0.85	A053420L	A053314F	A053814	A053914F	A049716F	
0.081	2.06	0.033	0.85	A053421L	A053315F	A053815	A053915F	A049717F	
0.090	2.29	0.033	0.85	A053422L	A053316F	A053816	A053916F	A049718F	
0.100	2.54	0.033	0.85	A053423L	A053317F	A053817	A053917F	A049719F	
0.110	2.79	0.033	0.85	A053424L	A053318F	A053818	A053918F	A049720F	
0.125	3.17	0.033	0.85	A053425L					



Saint-Gobain Performance Plastics Corporation
2664 Gilchrist Road, Akron, OH 44305
330-798-9240 • 800-798-1554 • 330-798-6968
www.tygon.com

Saint-Gobain Performance Plastics Corporation assumes no obligations or liability for any advice furnished by it, or for results obtained with respect to those products. All such advice is given and accepted at the buyer's risk.

©Saint-Gobain Performance Plastics Corporation